

ORDER INFORMATION

CODE : DL0901 - R1 - 1 X 8 ML + R2 - 1 X 2 ML
DL0902 - R1 - 1 X 20 ML + R2 - 1 X 5 ML


INTENDED USE :

This reagent kit is intended for "*in vitro*" quantitative determination of Creatine Kinase (CK-NAC) activity in serum based upon IFCC and DGKC recommendations.

CLINICAL SIGNIFICANCE :

Creatine kinase (CK) is an enzyme, which is found primarily in skeletal muscle, cardiac muscle and brain tissue. Elevated levels of CK are associated with myocardial infarction, various muscle disorders and diseases such as progressive Duchenne-type muscular dystrophy. In myocardial infarction, peak CK levels occur 24 to 36 hours after onset of chest pain and depending on the extent of damage can reach more than 10 times normal levels.

PRINCIPLE :



CK = Creatine kinase

HK = Hexokinase

G-6-P = Glucose-6-phosphate

G-6-PDH = Glucose-6-phosphate-dehydrogenase

REAGENT COMPOSITION :

Reagent 1 : Enzyme Reagent 1

Reagent 2 : Enzyme Reagent 2

MATERIALS REQUIRED BUT NOT PROVIDED :

- Clean & Dry Glassware.
- Micropipettes & Tips.
- Colorimeter or Bio-Chemistry Analyzer.

SAMPLES :

Serum free of hemolysis. Heparin or EDTA plasma.

WORKING REAGENT PREPARATION & STABILITY :

Mix 4 Volume of Enzyme Reagent 1, with 1 Volume of Reagent Enzyme 2. Working Reagent is stable for 30 days at 2 - 8°C.

GENERAL SYSTEM PARAMETERS :

Reaction type	Kinetic Reaction
Wave length	340 nm
Light Path	1 Cm
Reaction Temperature	37°C
Blank / Zero Setting	With Distilled Water
Reagent Volume	1ml
Sample Volume	50 µl
Lag / Delay Time	120 Sec.
Read Time	180 Sec.
Interval Time	60 Sec.
Factor	4127
Low Normal at 37°C	24 U/l
High Normal at 37°C	195 U/l
Linearity	1500 U/l
Max. Δ Abs / Min	0.363

ASSAY PROCEDURE :

Working Reagent	1000 ul
Sample	50 ul

Mix and after 120 second incubation, measure the change in absorbance every minute during 3 minutes at 37°C.

Determine the ΔA/min.

CALCULATION :

Creatine Kinase (CK-NAC) activity (U/l) = ΔA/min. x 4127

LINEARITY :

Reagent is Linear up to 1500 U/l

Dilute the sample appropriately and re-assay if Creatine Kinase (CK-NAC) activity exceeds 1500 U/l or ΔAbs / min Exceeds 0.363.

Multiply result with dilution factor.

REFERENCE NORMAL VALUE :

Female : 24-170 U/l

Male : 24-195 U/l

QUALITY CONTROL :

For accuracy it is necessary to run known controls with every assay.

LIMITATION & PRECAUTIONS :

1. Storage conditions as mentioned on the kit to be adhered.
2. Do not freeze or expose the reagents to higher temperature as it may affect the performance of the kit.
3. Before the assay bring all the reagents to room temperature.
4. Avoid contamination of the reagent during assay process.
5. Use clean glassware free from dust or debris.

BIBLIOGRAPHY :

Young et al., Clin. Chem., 21:10 (1975) Moren L.G., Clin. Chem., 23:1569 (1977)



DELTA LAB

1116, Bhadrakali Compound, Off Mumbai - Goa Highway, At Post Zarap,

Tal. : Kudal, Dist. Sindhudurg, Maharashtra - 416 510, INDIA

Tel / Fax : 0 23 62 - 23 20 30, Email : delta@deltalab.in, Web : www.deltalab.in